

NASOPHARYNGOLARYNGOSCOPE

Akino
Surgery Simplified



Nasopharyngolaryngoscope

Flexible Nasopharyngolaryngoscope

Product features

- HD optics in 3.4mm.
- Sturdy, long life optics as there is no working channel which reduces optics life.
- Inspection of nasal cavities, the pharynx and larynx upto vocal chords.
- Control lever with UP and DOWN direction for angular deflection of 130°
- Soakable, ETO Sterilizable
- Supplied with leakage tester
- Cold light cable connection
- Compatible for fiber optic cable of Akino or Storz type cable
- Eyepiece compatible with standard camera coupler

Technical Specifications

- Angle of view : 85°
- Outer diameter of the insertion sheath : 3.4mm
- Working length : 300mm
- Angular deflection : 130° up / down
- Protection class : IP1
- Medical product class : I
- Applied standard : IEC 60601-1
- Certification : CE



Operating Conditions -

- Temperature : 10°C to 55°C
- Relative humidity : 5 - 95%
- Air pressure : 50 - 106Kpa

Article No. - 2034.3000F

Accessories and Spare Parts -

- Pressure tester Art. No. : 2D00.NLPT
- ETO cap Art. No. : 2D00.NLET
- Wolf adapter Art. No. : 2D00.NLWA



Akino Medical Pte Ltd.

Japan

Level 9 Ariake Frontier Building Tower B,
3-7-26 Ariake, Koto-ku, Tokyo - 135-0063, Japan
Tel : +81 35530 8429. E-mail : info.jp@akinomedical.com

Singapore

362, Upper Paya Lebar Road, #04-09,
Da Jin Factory building, Singapore - 534963
Tel : + 65 68087810. E-mail : info.sg@akinomedical.com